

ID	Requirements	Related Use Case	Fulfilled By	Description	Implemented By	Tested By
1	A light on the machine and/or a symbol on the screen that indicates an active pulse reading	UC-01 UC-07	MainWindow	The moment the user starts a new session a light appears to signify that the machine is working according to adjustable settings	QTwidgets, readData()	readData()
2	The screen contains the menu options and the display graph. There are eight buttons: an off/on button for the device, a menu button, a standard back button which will return the user to the menu, four arrow buttons) and a selector in the centre which also functions as a start/stop button in session mode.	UC-01, UC-02 UC-03	MainWindow, Menu	All the buttons appear when the user turns on the heartwave.	QTframe, QPushButton, QTwidgets, readData(), navigateBack() mainMenu()	readData()
3	The device has a led light that changes to red, blue or green to indicate coherence level	UC-01, UC-07	MainWindow	At the top of Outer Frame, a widget changes colour based on the coherence level	QTwidgets, readData(), startSession(),	readData()
4	The menu options are displayed as default on the session screen. There is an option at the top to start a new session. The menu could consist of the following options: start new session, settings, log/history	UC-01, UC-03, UC-07	Menu, MainWindow	When the user powers on the device the main menu is the default screen	mainMenu(), initializeMenu(), updateMenu()	power_button()

ID	Requirements	Related Use Case	Fulfilled By	Description	Implemented By	Tested By
5	Session screen must display the main HRV graph (HR vs time) with key metric	UC-01 UC-03 UC-07	MainWindow	A Graph to display the coherence score	readData(), startSession(), addPoint()	startSession(), readData()
6	On the device there should be a light that changes to red, blue or green indicating low, medium or high coherence, depending on the challenge level	UC-07	MainWindow	At the top of Outer Frame, a widget changes colour based on the coherence, challenge level	readData()	startSession(), readData()
7	The metrics on the screen include the current coherence score (numerical value), length (duration of session), achievement (total sum of coherence scores sampled every 5 seconds	UC-01 UC-03	MainWindow	A new coherence, achievement score is displayed ever 5 seconds as well the length of session	readData(), updateValue()	startSession(), readData(), addPoint()
8	A breath pacer in the form of a strip of lights on the machine itself, or a ball going back and forth on the session screen ,default set at one breath every 10 seconds, adjustable in settings	UC-01 UC-04	MainWindow, Menu	Coloured widgets to display the coherence level picked by the user	navigateSubMenu(), navigateBack()	readData()
9	The settings tab includes challenge level (optional) and breath pacer setting	UC-01 UC-04 UC-05	Menu, MainWindow	There is an option to pick breath pacer level	navigateSubMenu(), navigateBack()	initializeMainMenu()
10	There are 4 challenge levels for coherence, from beginner to advanced, for the user to choose	UC-01 UC-05	Menu, MainWindow	There is an option to pick challenge	navigateSubMenu(), navigateBack()	initializeMainMenu()

ID	Requirements	Related Use Case	Fulfilled By	Description	Implemented By	Tested By
11	The breath pacer, 1-30 seconds, increases time interval between each breath, default at 10 seconds	UC-01 UC-04	Menu, MainWindow	Pick a breath pacer level from settings	navigateSubMenu(), navigateBack()	initializeMainMenu()
12	When the user ends a session a summary view will appear that includes the following information: challenge level (optional), percentage of time in different coherence levels (low, medium and high), average coherence, length of session, achievement score, entire HRV graph	UC-01 UC-03 UC-06 UC-08	MainWindow, qcustomplot	When a session ends, save all the data to History	endSession(), updateHistory()	endSession(), navigateSubMenu(),
13	The menu contains a log or history tab of all sessions, with dates, when selected show the summary view, as well as the ability to delete a session	UC-01 UC-08	Menu, MainWindow	Sessions are available to be chosen to review	mainMenu(), updateHistory(), clearHistory()	updateHistory(), navigateSubMenu(),
14	An option to reset, wipe all data and restore the device to the initial install condition	UC-01 UC-08	Menu, MainWindow	In the History menu, there is an option to clear/reset device	clearHistory()	updateHistory(), clearHistory()
15	There is a battery charge indicator on the session screen	UC-01 , UC-02	MainWindow	There is a widget displaying the battery life of the HeartWave	updateValues(), chargeBattery()	chargeBattery()
16	A beep goes off when a new coherence level is reached	UC-01 UC-06	MainWindow	Notifies the user when the coherence level changes	readData()	readData(), updateValues()